

## Protocol: DNA Purification from Yeast Using the Gentra Puregene Yeast/Bact. Kit

This protocol is for purification of genomic DNA from fresh or frozen samples of 1 ml overnight yeast cultures (approximately  $1-2 \times 10^8$  cells) using the Gentra Puregene Yeast/Bact. Kit.

### Things to do before starting

- Preheat water baths to 37°C for use in steps 6 and 20 and 65°C for use in step 21 of the procedure.
- Frozen yeast samples should be thawed and equilibrated to room temperature (15–25°C) before beginning the procedure.

### Procedure

1. Prepare an overnight culture containing  $1-2 \times 10^8$  cells.
2. Transfer 1 ml of the cell suspension to a 1.5 ml microcentrifuge tube on ice.
3. Centrifuge for 5 s at 13,000–16,000  $\times g$  to pellet cells.
4. Carefully discard the supernatant by pipetting or pouring.
5. Add 300  $\mu$ l Cell Suspension Solution, and pipet up and down.
6. Add 1.5  $\mu$ l Lytic Enzyme Solution, and mix by inverting 25 times. Incubate for 30 min at 37°C.
7. Centrifuge for 1 min at 13,000–16,000  $\times g$  to pellet cells.
8. Carefully discard the supernatant by pipetting or pouring.
9. Add 300  $\mu$ l Cell Lysis Solution, and pipet up and down to lyse the cells.
10. Add 100  $\mu$ l Protein Precipitation Solution, and vortex vigorously for 20 s at high speed.
11. Centrifuge for 3 min at 13,000–16,000  $\times g$ .  
The precipitated proteins should form a tight pellet. If the protein pellet is not tight, incubate on ice for 5 min and repeat the centrifugation.
12. Pipet 300  $\mu$ l isopropanol into a clean 1.5 ml microcentrifuge tube and add the supernatant from the previous step by pouring carefully.  
Be sure the protein pellet is not dislodged during pouring.
13. Mix by inverting gently 50 times.
14. Centrifuge for 1 min at 13,000–16,000  $\times g$ .  
The DNA may be visible as a small white pellet.
15. Carefully discard the supernatant, and drain the tube by inverting on a clean piece of absorbent paper, taking care that the pellet remains in the tube.

16. Add 300  $\mu$ l of 70% ethanol and invert several times to wash the DNA pellet.
17. Centrifuge for 1 min at 13,000–16,000  $\times g$ .
18. Carefully discard the supernatant. Drain the tube on a clean piece of absorbent paper, taking care that the pellet remains in the tube. Allow to air dry for 5 min.  
The pellet might be loose and easily dislodged. Avoid over-drying the DNA pellet, as the DNA will be difficult to dissolve.
19. Add 100  $\mu$ l DNA Hydration Solution and vortex for 5 s at medium speed to mix.
20. Add 1.5  $\mu$ l RNase A Solution, and mix by vortexing by 1 s. Pulse spin to collect liquid, and incubate at 37°C for 15–60 min.
21. Incubate at 65°C for 1 h to dissolve the DNA.
22. Incubate at room temperature (15–25°C) overnight with gentle shaking. Ensure tube cap is tightly closed to avoid leakage. Samples can then be centrifuged briefly and transferred to a storage tube.